

CEN FISH probe

Volume: 200ul

	p-arm	q-arm
1	CEN1p/780kb	CEN1q/600kb
2	CEN2p/670kb	CEN2q/580kb
3	—	CEN3q/500kb
4	CEN4p/630kb	CEN4q/650kb
5	CEN5p/490kb	CEN5q/700kb
6	CEN6p/660kb	—
7	CEN7p/700kb	CEN7q/800kb
8	CEN8p/520kb	CEN8q/520kb
9	—	CEN9q/470kb
10	CEN10p/640kb	CEN10q/450kb
11	CEN11p/630kb	CEN11q/680kb
12	CEN12p/530kb	CEN12q/580kb
13	—	CEN13q/550kb
14	—	CEN14q/400kb
15	—	CEN15q/680kb
16	CEN16p/360kb	CEN16q/700kb
17	—	CEN17q/540kb
18	CEN18p/600kb	CEN18q/660kb
19	CEN19p/570kb	CEN19q/430kb
20	CEN20p/400kb	CEN20q/680kb
21	—	CEN21q/520kb
22	—	CEN22q/620kb
X	CENXp/550kb	
Y	—	CENYq/470kb